




Miltenyi Biotec

MACS® Technology

Magnetic mitochondria isolation



Functional mitochondria

High purity

Outstanding yield

Two-hour workflow



► miltenyibiotec.com/mitochondria

Mitochondria Isolation Kit, human

For 25 isolations from human tissue or cells




Figure 1: Principle of MACS Technology for mitochondria isolation

The **Mitochondria Isolation Kit, human** enables isolation of intact, viable mitochondria from human cells or tissue in less than two hours. The kit's protocol is based on the renowned MACS® Technology that combines high yield with purity and integrity of the isolated mitochondria (figs. 2 and 3). After cell lysis, the translocase of the outer mitochondrial membrane 22 (TOM22) is targeted by Anti-TOM22 MicroBeads, human to enable magnetic isolation of functional human mitochondria. See figure 1 for details.

Functional mitochondria preparation




Figure 2: Mitochondria were prepared using the Mitochondria Isolation Kit, human (MACS) or by differential centrifugation (DC) and were then incubated for 1 h with ³⁵S-labeled TFAM. Import was studied by SDS-PAGE and fluoro-graphy. Lane 1: Control lysate of *in vitro* translated TFAM. Lane 2-7: Proteins from sedimented mitochondria. P= Precursor protein. M= Imported mature protein. The data suggest that the function of mitochondria preparations isolated with the Mitochondria Isolation Kit, human is fully maintained. (Courtesy of Dr. H.-T. Hornig-Do, Newcastle, UK)

Purest mitochondria preparation




Figure 3: Mitochondria were prepared by different protocols. Amounts of COX I (mitochondria) and KDEL (endoplasmic reticulum) in the mitochondrial fractions obtained using differential centrifugation (DC), Mitochondria Isolation Kit, human (MACS), and ultracentrifugation (UC) were analyzed by Western blotting. Compared to DC and UC, mitochondria isolated with the Mitochondria Isolation Kit, human have no detectable contamination with ER. (Courtesy of Dr. H.-T. Hornig-Do, Newcastle, UK)

Mitochondria Isolation Kit, mouse tissue

For 25 isolations from mouse tissue




Figure 4: Integrated tissue dissociation followed by mitochondria isolation.

The **Mitochondria Isolation Kit, mouse tissue** is intended for isolation of mitochondria from any mouse tissue such as brain, liver, or skeletal muscle. Tissues can be dissociated with virtually any protocol although the use of the gentleMACS™ Dissociator is preferred because this will save time and ensures standardization of tissue dissociation and reproducible results (see fig. 4). Similar to the Mitochondria Isolation Kit, human, this kit allows for enrichment of functional and pure mitochondria (figs. 5 and 6).

Enrichment of functional mitochondria




Figure 5: Mitochondria were prepared from various mouse tissues using the Mitochondria Isolation Kit, mouse tissue. The mitochondria preparations from muscle, liver and brain as well as the flow-through fractions were analyzed by a coupled enzymatic citrate synthase assay. The data indicate that citrate synthase activity, normalized to the activity per mg protein assayed in the flow-through fraction, is highly enriched in the eluate.

High purity liver mitochondria





Figure 6: After isolation from liver using the Mitochondria Isolation Kit, mouse tissue, equal amounts of mitochondrial protein lysates were analyzed by Western blotting against TOM20. PEX1, a cytoplasmic protein anchored to peroxisomal membranes, served as a control. The results show an effective enrichment of mitochondria and no contamination of the mitochondria fraction by peroxisomes.

Product overview

Product	Capacity	Components	Order no.
Mitochondria isolation from human tissues or cells			
Mitochondria Isolation Kit, human	25 isolations	1.25 mL Anti-TOM22 MicroBeads, human, buffers, LS Columns	130-094-532
Mitochondria MidiMACS Starting Kit, human	25 isolations	1 Mitochondria Isolation Kit, 1 MidiMACS™ Separator, 1 MACS MultiStand	130-094-872
Mitochondria QuadroMACS Starting Kit, human	25 isolations	1 Mitochondria Isolation Kit, 1 QuadroMACS™ Separator, 1 MACS MultiStand	130-094-833
Mitochondria isolation from mouse tissue			
Mitochondria Isolation Kit, mouse tissue	25 isolations	1.25 mL Anti-TOM22 MicroBeads, mouse, buffers, LS Columns, Pre-Separation Filters, 30 µm	130-096-946
Mitochondria MidiMACS Starting Kit, mouse tissue	25 isolations	Mitochondria Isolation Kit, MidiMACS Separator, 1 MACS MultiStand	130-097-039
Mitochondria QuadroMACS Starting Kit, mouse tissue	25 isolations	Mitochondria Isolation Kit, QuadroMACS Separator, 1 MACS MultiStand	130-097-040
Related products			
gentleMACS Dissociator	1–2 samples in parallel	Instrument, 25 C Tubes, 25 M Tubes, power cord, user manual	130-093-235
gentleMACS Octo Dissociator	up to 8 samples in parallel	Instrument, power cord, user manual	130-095-937
Neural Tissue Dissociation Kit (P)	50 digestions	2.5 mL Solution 1, 2x50 mL Solution 2, 1.5 mL Solution 3, 1 vial Solution 4, 1 mL Storage Buffer	130-092-628
Neural Tissue Dissociation Kit (T)	50 digestions	10 mL Solution 1, 2x50 mL Solution 2, 1.5 mL Solution 3, 1 vial Solution 4, 1 mL Storage Buffer	130-093-231
Tumor Dissociation Kit, human	25 digestions	2 vials Solution 1, 1 vial Solution 2, 1 vial Solution 3, 1 mL Reconstitution Buffer	130-095-929
Tumor Dissociation Kit, mouse	50 digestions	2 vials Solution 1, 1 vial Solution 2, 1 vial Solution 3, 1 mL Reconstitution Buffer	130-096-730
Brain Tumor Dissociation Kit (P), human	25 digestions	1.25 mL Solution 1, 2x50 mL Solution 2, 1.5 mL Solution 3, 1 vial Solution 4, 1 mL Storage Buffer	130-095-942
Brain Tumor Dissociation Kit (T), human	25 digestions	4 mL Solution 1, 2x50 mL Solution 2, 1.5 mL Solution 3, 1 vial Solution 4, 1 mL Storage Buffer	130-095-939



Testimonials and References

“Well-coupled mitochondria are a prerequisite to achieving reliable reproducible results in nearly all functional assays...mitochondria isolated with the MACS protocol were found to be better coupled than those obtained using DC...”

Hornig-Do, H.T. *et al.* (2009) *Anal. Biochem.* 389: 1–5.

“In order to maximize mitochondria enrichment efficiency and minimize contaminations from other organelles, we employed a newly developed method for mitochondria isolation based on superparamagnetic microbeads conjugated to anti-TOM22 antibody. The protocol is fast, reproducible, and standardized, resulting in mitochondria of high purity, with minimal contamination from cytoskeleton, cytosol, Golgi apparatus, endosome, endoplasmic reticulum, and nucleus.”

Guo, T. *et al.* (2010) *Mol. Cell. Proteomics* 9: 2629–2641.

References

Hornig-Do, H.T. *et al.* (2009) Isolation of functional pure mitochondria by superparamagnetic microbeads. *Anal. Biochem.* 389: 1–5.

Minet, A.D. and Gaster, M. (2011) The dynamic equilibrium between ATP synthesis and ATP consumption is lower in isolated mitochondria from myotubes established from type 2 diabetic subjects compared to lean control. *Biochem. Biophys. Res. Comm.* 409: 591–595.

Minet, A.D. and Gaster, M. (2010) ATP synthesis is impaired in isolated mitochondria from myotubes established from type 2 diabetic subjects. *Biochem. Biophys. Res. Comm.* 402: 70–74.

Guo, T. *et al.* (2010) Quantitative proteomics discloses MET expression in mitochondria as a direct target of MET kinase inhibitor in cancer cells. *Mol. Cell. Proteomics* 9: 2629–2641.

Eriksen, M.B. *et al.* (2011) Intact Primary Mitochondrial Function in Myotubes Established from Women with PCOS. *J. Clin. Endocrinol. Metab.* 96: E1298–1302.

Barrey, E. *et al.* (2011) Pre-microRNA and mature microRNA in human mitochondria. *PLoS ONE* 6: e20220.

Bandiera, S. *et al.* (2011) Nuclear outsourcing of RNA interference components to human mitochondria. *PLoS ONE* 6: e20746.

Sacchi, S. *et al.* (2011) Evidence for the interaction of d-amino acid oxidase with pLG72 in a glial cell line. *Mol. Cell. Neurosc.* 48: 20–28.



Miltenyi Biotec

Germany/Austria/ Switzerland

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

USA/Canada

Miltenyi Biotec Inc.
2303 Lindbergh Street
Auburn, CA 95602, USA
Phone 800 FOR MACS
Phone +1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

Australia

Miltenyi Biotec
Australia Pty. Ltd.
Unit 16A, 2 Eden Park Drive
North Ryde, NSW 2113,
Australia
Phone +61 2 8877 7400
Fax +61 2 9889 5044
macs@miltenyibiotec.com.au

Benelux

Miltenyi Biotec B.V.
Schipholweg 68 H, 2316 Leiden
The Netherlands
macs@miltenyibiotec.nl

Customer service Netherlands

Phone 0800 4020120
Fax 0800 4020100

Customer service Belgium

Phone 0800 94016
Fax 0800 99626

Customer service Luxembourg

Phone 800 24971
Fax 800 24984

China

Miltenyi Biotec GmbH
Shanghai Office
Rm. 2309-2310,
No. 319 Xianxia Rd.
Shanghai 200051, P.R. China
Phone +86 21 62351005
Fax +86 21 62350953
macs@miltenyibiotec.com.cn

France

Miltenyi Biotec SAS
10 rue Mercœur
75011 Paris, France
Phone +33 1 56 98 16 16
Fax +33 1 56 98 16 17
macs@miltenyibiotec.fr

Italy

Miltenyi Biotec S.r.l.
Via Persicetana, 2/D
40012 Calderara di Reno (BO)
Italy
Phone +39 051 6 460 411
Fax +39 051 6 460 499
macs@miltenyibiotec.it

Japan

Miltenyi Biotec K.K.
Nittsu-Eitai Building 5F
16-10 Fuyuki, Koto-ku,
Tokyo 135-0041, Japan
Phone +81 3 5646 8910
Fax +81 3 5646 8911
macs@miltenyibiotec.jp

Singapore

Miltenyi Biotec
Asia Pacific Pte Ltd.
100 Beach Road
#28-06 to 28-08 Shaw Tower
Singapore 189702
Phone +65 6238 8183
Fax +65 6238 0302
macs@miltenyibiotec.com.sg

Spain

Miltenyi Biotec S.L.
C/Luis Buñuel 2
Ciudad de la Imagen,
28223 Pozuelo de Alarcón
(Madrid), Spain
Phone +34 91 512 12 90
Fax +34 91 512 12 91
macs@miltenyibiotec.es

United Kingdom

Miltenyi Biotec Ltd.
Almac House, Church Lane
Bisley, Surrey GU24 9DR, UK
Phone +44 1483 799 800
Fax +44 1483 799 811
macs@miltenyibiotec.co.uk

www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, gentleMACS, MidiMACS, QuadroMACS, and the MACS Logo are registered trademarks of Miltenyi Biotec GmbH. Copyright © 2011 Miltenyi Biotec GmbH. All rights reserved.