Miltenyi Biotec

# MACSmix™ Tube Rotator

## User manual

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark and gentleMACS, MACSiBead, and MACSmix are trademarks of Miltenyi Biotec GmbH.

Copyright © 2009 Miltenyi Biotec GmbH. All rights reserved.

No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec GmbH; however, notwithstanding the foregoing, the owners of the MACSmix™ Tube Rotator may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization.

Maximal care has been taken by Miltenyi Biotec GmbH in the preparation of this user manual. However, Miltenyi Biotec GmbH shall not be liable for any technical or editorial errors or omissions contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided "as is" without warranty of any kind, and is subject to change without notice.

# **MACSmix<sup>TM</sup>**

# **Tube Rotator**

## **User manual**

Miltenyi Biotec GmbH  
Friedrich-Ebert-Straße 68  
D – 51429 Bergisch Gladbach  
Germany

Phone +49 2204 8306-0  
Fax +49 2204 85197  
macstec@miltenyibiotec.de

[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

**Version 06**

# Contents

	<b>Warnings and precautions</b>	<b>6</b>
	<b>Warnhinweise und Vorsichtsmaßnahmen</b>	<b>8</b>
	<b>Avertissements et précautions</b>	<b>10</b>
	<b>Avvertimenti e precauzioni</b>	<b>12</b>
	<b>Advertencias y precauciones</b>	<b>14</b>
	<b>Technical service</b>	<b>16</b>
<b>Chapter 1</b>	<b>Introduction</b>	<b>17</b>
	1.1 Description of the MACSmix™ Tube Rotator	17
	1.2 Applications	17
	1.3 Supplied components	17
<b>Chapter 2</b>	<b>Installation</b>	<b>18</b>
<b>Chapter 3</b>	<b>Instructions for use</b>	<b>19</b>
	3.1 Running modes	19
	3.2 The racks	20
	3.3 Function of the keys	22
	3.4 Operating the MACSmix™ Tube Rotator	23
<b>Chapter 4</b>	<b>Maintenance</b>	<b>24</b>
	4.1 Charging the batteries	24
	4.2 Changing the batteries	25
	4.3 Cleaning the MACSmix™ Tube Rotator	26
	4.4 Decontamination of the MACSmix™ Tube Rotator	27
<b>Chapter 5</b>	<b>Troubleshooting</b>	<b>28</b>
<b>Chapter 6</b>	<b>Technical data</b>	<b>29</b>

# Warnings and precautions

## **CAUTION**

Documentation needs to be consulted before proceeding with installation and operation of the system.

## **WARNING**

A hazard exists if the operating instructions are not followed.

## **CAUTION**

The MACSmix™ Tube Rotator is for research applications only.

## **CAUTION**

The MACSmix Tube Rotator should only be used in the range of temperatures between 2 °C and 42 °C. The device may be damaged if it is used at higher or lower temperatures.

## **CAUTION**

The MACSmix Tube Rotator is intended for indoor use only. Do not use the instrument in locations classified as hazardous or with hazardous materials for which the equipment is not designed.

## **CAUTION**

The safety features of the MACSmix Tube Rotator may be impaired if the equipment is used with accessories not provided or recommended by Miltenyi Biotec, or used in a manner not specified by Miltenyi Biotec.

## **CAUTION**

A hazard may occur if the MACSmix Tube Rotator is used to mix materials classified as inflammable or where the transfer of mechanical energy to glass apparatus could lead to breakage.

## **WARNING**

As electric devices can cause an electric shock, special care must be taken while handling fluids. Please ensure that fluids do not enter the interior part of the device. Unplug the power cord before cleaning the MACSmix Tube Rotator.

**CAUTION**

Check the line voltage on the power supply unit, it must correspond to the voltage indicated. The electrical components of the MACSmix Tube Rotator may be damaged if the power supply and line voltages are not compatible.

**CAUTION**

Do not fit a battery of the wrong type in the battery compartment. An explosion or fire hazard could occur. Do not fit non-rechargeable cells in the battery compartment.

**NOTICE**

Use only NiMH rechargeable batteries type AA with 2000–4000 mAh. Check the correct polarity, otherwise the MACSmix Tube Rotator can be damaged.

**NOTE:** NiMH rechargeable batteries have to be disposed of as hazardous waste.

**NOTE:** Electromagnetic Compatibility (EMC): Changes or modifications of the equipment unless expressly approved by Miltenyi Biotec may void your authority to operate the equipment pursuant to Part 15 of the Federal Communications Commission (FCC) Rules.

## Limited warranty

Should the MACSmix Tube Rotator be used in a manner not specifically described in this manual, all warranties will be null and void.

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# Technical service

**For technical questions, please contact your local distributor  
or our technical support team:**

Miltenyi Biotec GmbH  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany

Phone +49 2204 8306 830  
Fax +49 2204 8306 89

[macstec@miltenyibiotec.de](mailto:macstec@miltenyibiotec.de)  
[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

## 1.1 Description of the MACSmix™ Tube Rotator

The MACSmix™ Tube Rotator is a versatile tube rotator that runs on rechargeable batteries, working independently from permanent power supply. It can be used at a temperature range of 2 °C to 42 °C and can be placed in an incubator as well as in a refrigerator. The MACSmix Tube Rotator can be equipped with tubes from 0.5 mL to 50 mL in size. The device operates at three different speeds or at set intervals. Thus, the MACSmix Tube Rotator is an ideal tube rotator for a variety of applications.

## 1.2 Applications

### For example:

- Secretion period during MACS® Cytokine Secretion Assays
- Enzymatic steps during tissue dissociation
- Short time culturing of cells in tubes in an incubator
- Loading of MACSiBead™ Particles with biotinylated antibodies for cell activation
- Magnetic labeling of cells for depletion using MACSiBead Particles
- Resuspending blood or other cell suspensions
- Immunoprecipitations
- Binding of molecules, e.g., recombinant proteins to an affinity matrix
- Extractions and phase separations
- Dissolving substances

### CAUTION

If normal use involves the handling of hazardous substances, instruction should be given on correct use and safety provisions.

## 1.3 Supplied components

- **1 MACSmix™ Tube Rotator (# 130-090-753)**
- 4 rechargeable NiMH batteries, type AA
- 1 power supply for 100–240 V, 50–60 Hz
- 1 small rack for tubes from 0.5 mL to 2 mL in size
- 1 large rack for tubes from 5 mL to 50 mL in size
- 1 screwdriver
- 1 User manual

# 2 Installation

**Notice:** Insert the plug specific for your country into the power supply socket.

1. Open the lid on the base of the device by using the screwdriver. Insert the 4 batteries into the battery holder in the indicated direction. Close the lid by using the screwdriver.

**Note:** Wrong orientation of the batteries can damage the device

**Note:** The capacity selector is already set to the correct position for the enclosed batteries. For instructions on changing the batteries refer to section 4.2 Changing the batteries.

2. Connect the power cable into the outlet and plug the power supply into the socket of the MACSmix™ Tube Rotator.

**Note:** Using other than the original power supply may damage the device.

3. Charge the batteries for 2.5–4 hours depending on the capacity of the batteries, during this time the green LED above the OFF-key will flash.

Battery capacity	Charging time
2000 mAh	2.5 h
3000 mAh	3 h
4000 mAh	4 h

Once the recharging process has finished, the green light stops flashing. The batteries cannot be overcharged. The MACSmix Tube Rotator may be used while charging.

**Note:** Make sure that the batteries are fully charged before using the MACSmix Tube Rotator for the first time.

## 3.1 Running modes

### **Permanent run on the power supply**

The MACSmix Tube Rotator can be operated permanent with the power supply. Overcharging the batteries is prevented by a microprocessor.


### **Run on rechargeable batteries**

The MACSmix Tube Rotator can be used without a permanent connected power supply. The device can be put into a refrigerator or an incubator at temperatures between 2 °C and 42 °C for operation. Fully recharged batteries with minimum of 2000 mAh capacity enable the device to operate for at least 48 hours at room temperature. This time could be shorter at higher or lower temperatures.

## 3.2 The racks


### Small rack for 0.5 mL to 2 mL tubes

- The small rack can be equipped with up to 16 tubes of 1.5 mL to 2 mL in size on one side, or with up to 16 tubes of 0.5 mL in size by using the other side of the rack.
- The small rack is then placed in the MACSmix™ Tube Rotator fitting the coupling of the rack (1) into the bearing parts of the MACSmix Tube Rotator (2).

### Large rack for 5 mL to 50 mL tubes

- The large rack can be used for 8 tubes of 5 mL up to 15 mL in size, or for 4 tubes of 50 mL in size.
- The large rack can be adjusted to the appropriate tube length of the tubes by unscrewing the movable spring loaded disc **(2)** and securing it in the new position by tightening the screw **(3)**.
- The tubes are fixed with the spring loaded disc. The rack should be loaded symmetrically. For sufficient mixing, the tubes need to be placed in a sloping position into the rack **(4)**.
- The large rack is then placed in the MACSmix Tube Rotator fitting the coupling of the rack **(1)** into the bearing parts of the MACSmix Tube Rotator.

### 3.3 Function of the keys

Press this key to switch off the MACSmix™ Tube Rotator. The LED flashes green while the batteries are being recharged and it stops flashing when the recharging process has finished. The LED flashes red when the remaining time of the batteries is less than 6 hours.

Press this key to operate the MACSmix Tube Rotator on permanent run at a speed of approximately 4 rpm. The green LED indicates this program.

Press this key to operate the MACSmix Tube Rotator on permanent run at a speed of approximately 8 rpm. The green LED indicates this program.

Press this key to operate the MACSmix Tube Rotator on permanent run at a speed of approximately 12 rpm. The green LED indicates this program.

Press this key to operate the MACSmix Tube Rotator at intervals of 30 seconds run at a speed of approximately 12 rpm and 1 minute pause. The green LED indicates this program.

Press this key to operate the MACSmix Tube Rotator at intervals of 30 seconds run at a speed of approximately 12 rpm and 5 minutes pause. The green LED indicates this program.

## 3.4 Operating the MACSmix™ Tube Rotator

1. For operation without permanent connected power supply, make sure that the batteries are charged (see 4.1).

For operation with permanent connected power supply, connect the power cable into the outlet and plug the power supply into the socket of the MACSmix™ Tube Rotator.

Make sure that the power cable is easily accessible and located as close to the operator of the instrument as possible.

**Note:** If it is necessary to disconnect the power supply unplug the power cable from the power outlet.

2. Choose the suitable rack according to the size of the tubes and place the tubes into the rack (see 3.2).
3. Install the rack in the MACSmix Tube Rotator by placing the coupling of the rack into the bearing parts of the MACSmix Tube Rotator (see 3.2).
4. Choose the appropriate program by pressing the respective key (see 3.3).
5. Stop the operation of the MACSmix Tube Rotator by pressing the OFF-key (see 3.3).

# 4 Maintenance

## 4.1 Charging the batteries

1. The batteries of the MACSmix™ Tube Rotator should only be recharged with the enclosed power supply.

**Note:** Using other than the original power supply may damage the device.

2. Plug the main adapter of the transformer into the socket of the MACSmix Tube Rotator. Insert the transformer into the socket.
3. Charge the batteries for 2.5–4 hours depending on the capacity of the batteries, during this time the green LED above the OFF-key will flash.

Battery capacity	Charging time
2000 mAh	2.5 h
3000 mAh	3 h
4000 mAh	4 h


Once the recharging process has finished, the green light stops flashing.

**Note:** The batteries cannot be overcharged. The MACSmix Tube Rotator may be used while recharging.

## 4.2 Changing the batteries

1. For replacement of batteries that are no longer rechargeable disconnect the device from the power supply. Open the lid on the base of the device by using the screwdriver.
2. Remove the old batteries from the battery holder. Set the capacity selector switch to the correct position. The capacity is printed on the batteries.

**Note:** The capacity selector is already set to the correct position for the enclosed batteries.

Battery capacity	Capacity selector
2000–2600 mAh	Position 1
2600–4000 mAh	Position 2

Insert the 4 new batteries into the battery holder as indicated. Close the lid by using the screwdriver.

**Note:** Wrong orientation of the batteries can damage the device.

**Note:** NiMH batteries have to be disposed of as hazardous waste.

### 4.3 Cleaning the MACSmix™ Tube Rotator

Switch off and unplug the MACSmix™ Tube Rotator before cleaning. Use only commercially available cleaning agents such as mild detergent or alcohol (70% ethanol) to clean the MACSmix device and the racks.

**Note:** Never use acetone!

## 4.4 Decontamination of the MACSmix™ Tube Rotator

Switch off and unplug the MACSmix™ Tube Rotator before decontamination.

- If the MACSmix Tube Rotator has been used to process biohazardous samples, it is recommended to decontaminate the instrument.
- If sample tubes have been in contact with biohazardous material, they should be autoclaved after use.
- In case of spillage, it is recommended to use a disinfectant that is appropriate for the potential pathogen, e.g., bleach, isopropyl alcohol, or 70% ethanol, to decontaminate surfaces with tissue or swabs.

**Note:** Dispose tissues and swabs appropriately. It is recommended to wear protective gloves, protective clothing, and safety glasses to prevent contact with skin and eyes.

# 5 Troubleshooting

<b>Error</b>	<b>Source of error</b>	<b>User action</b>	<b>See</b>
Racks do not move, but functional LED is on.	Bearing parts or coupling of the rack are soiled.	Clean the bearing parts and the couplings of the racks.	4.3
Racks do not move and functional LED does not work.	Batteries are empty.	Charge batteries.	4.1
OFF-key LED does not work while the batteries are charging.	Transformer has not been plugged in correctly.	Plug in transformer.	4.1
Running time is less than 48 hours, after fully charging the batteries.	Batteries are expired.	Change batteries.	4.2

<b>Dimensions</b>	W 281×D 96×H 85 (111 with large rack) mm
<b>Weight</b>	MACSmix Tube Rotator: 430 g Small rack: 220 g Large rack: 280 g
<b>Running modes</b>	Permanent mixing: 4 rpm, 8 rpm, or 12 rpm Interval mixing: 30 seconds running at 12 rpm and 1 minute or 5 minutes pause
<b>Tube capacity</b>	16 tubes with 0.5 mL, 1.5 mL, or 2 mL in size 8 tubes with 5 mL to 15 mL in size 4 tubes with 50 mL in size, or 4 gentleMACS™ Tubes
<b>Power supply</b>	Input voltage 100–240 VAC Frequency 50–60 Hz Power consumption 25 W Output voltage 12 VDC
<b>Electrical data</b>	Input voltage 12 VDC
<b>MACSmix Tube Rotator</b>	Power consumption 20 W
<b>Batteries</b>	4 NiMH rechargeable batteries 2000–4000 mAh type AA
<b>Standards</b>	The MACSmix Tube Rotator has been tested to and satisfies the requirements for EN 61010-1/A2:1998. Satisfies low voltage directive (LVD) specifications 73/23/EEC.
<b>Conditions for operation</b>	Temperature 2 °C to 42 °C Humidity 0% to 85% rel. Altitude maximal 2000 m
<b>Warranty</b>	12 months
<b>Protection class</b>	III
<b>Voltage</b>	Supply voltage fluctuations up to ±10% of the nominal voltage. Transient overvoltages present on the mains supply: Category II. The MACSmix Tube Rotator is suitable for rated pollution degree 2.

The MACSmix Tube Rotator meets the intent of the directive 2004/108/EC (electromagnetic compatibility) and 2006/95/EC (low voltage equipment). Compliance was demonstrated by conformance to the following standards which have been listed in the Official Journal of the European Communities:

**EMC:** EN 61326-1  
EN 61000-3-2  
EN 61000-3-3  
**Safety:** EN 61010-1  
EN 61010-2-051

**Compliance was demonstrated by conformance to the following FCC rules:**

FCC 47 CFR Ch.1 part 15, class B.

**Design and specification are subject to change without notice.**

**Order no. 130-090-753**

**For technical questions, please contact your local distributor or our technical support team in Germany:**

**macstec@miltenyibiotec.de**  
**Phone +49 2204 8306-830**  
**Fax +49 2204 8306-89**

**Miltenyi Biotec**

**Miltenyi Biotec GmbH**  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
Fax (orders) +49 2204 85197  
macs@miltenyibiotec.de

**Miltenyi Biotec Inc.**  
2303 Lindbergh Street  
Auburn, CA 95602  
USA  
Phone 800 FOR MACS,  
+1 530 888 8871  
Fax +1 530 888 8925  
macs@miltenyibiotec.com

**Miltenyi Biotec B.V. (Benelux)**  
macs@miltenyibiotec.nl  
**Customer service Netherlands**  
Phone 0800 4020120  
**Customer service Belgium**  
Phone 0800 94016  
**Customer service Luxembourg**  
Phone 800 24971

**Miltenyi Biotec SAS (France)**  
Phone +33 1 56 98 16 16  
macs@miltenyibiotec.fr  
**Miltenyi Biotec S.r.l. (Italy)**  
Phone +39 051 646 0411  
macs@miltenyibiotec.it

**Miltenyi Biotec Asia Pacific Pte. Ltd. (Singapore)**  
Phone +65 6238 8183  
macs@miltenyibiotec.com.sg

**Miltenyi Biotec S.L. (Spain)**  
Phone +34 91 512 12 90  
macs@miltenyibiotec.es

**Miltenyi Biotec K.K. (Japan)**  
Phone +81 3 5646 8910  
macs@miltenyibiotec.jp

**Miltenyi Biotec Ltd. (UK)**  
Phone +44 1483 799 800  
macs@miltenyibiotec.co.uk

**www.miltenyibiotec.com**

**Miltenyi Biotec Australia Pty. Ltd.**  
Phone +61 2 8877 7400  
macs@miltenyibiotec.com.au

**Miltenyi Biotec Trading (Shanghai) Co., Ltd. (China)**  
Phone +86 21 6235 1005  
macs@miltenyibiotec.com.cn

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use. MACS is a registered trademark and gentleMACS, MACSiBead, and MACSmix are trademarks of Miltenyi Biotec GmbH. Copyright © 2009 Miltenyi Biotec GmbH. All rights reserved.

130-091-017.06