



Human FGF-2 premium grade

10 µg	130-093-839	200 µg	130-093-841
50 µg	130-093-840	1000 µg	130-093-842
100 µg	130-093-564	2000 µg	130-093-843

Contents

1. Description

1.1 Background information

1.2 Applications

2. Reference

1. Description

Components	Human FGF-2, premium grade: Purified recombinant human fibroblast growth factor 2.
Sizes	10 µg, 50 µg, 100 µg, 200 µg, 1000 µg, 2000 µg.
Biological activity	The ED ₅₀ is <1.25 ng/mL* corresponding to a specific activity of >8×10 ⁵ I.U./mg.
Molecular mass	17 kDa.
Source	Produced in <i>E. coli</i> .
Product format	Lyophilized from a 0.2 µm filtered buffer solution.
Stabilizer	None (or trehalose and mannitol, 100 µg size).
Purity	>97% as determined by SDS-PAGE analysis.
Endotoxin level	Low endotoxin (<1 EU/µg cytokine) as determined by Limulus Amebocyte Lysate (LAL) assay.
Storage	Lyophilized Human FGF-2, premium grade should be stored at -20 °C. The expiration date is indicated on the vial label. Upon reconstitution aliquots should be stored at -20 °C. Avoid repeated freeze-thaw cycles.
Reconstitution	It is recommended to reconstitute lyophilized Human FGF-2 with deionized sterile-filtered water up to a final concentration of 100 µg/mL. ▲ Note: Addition of carrier protein, such as 0.1% bovine serum albumin (BSA) or human serum albumin (HSA) may have stabilizing effects. Further dilutions should be prepared with 1% BSA or HSA in phosphate-buffered saline (PBS).

* The ED₅₀ is determined by proliferation assay using 3T3 cells according to Robinson and Gaines-Das.¹ The proliferation assay was calibrated with the first international reference standard for human FGF-2 (NIBSC code 90/712) provided by the WHO/ National Institute for Biological Standards and Control.

1.1 Background information

Fibroblast growth factor 2 (FGF-2 or FGF2), also termed fibroblast growth factor basic (FGF-b) or basic FGF (bFGF), belongs to the family of heparin-binding growth factors. It functions as a wide-spectrum mitogenic, angiogenic, and neurotrophic factor and stimulates the proliferation of a wide variety of cells including mesenchymal, neuroectodermal, and endothelial cells. FGF-2 has been implicated in a multitude of physiological and pathological processes, including limb development, angiogenesis, wound healing, and tumor growth.

1.2 Applications

- Human FGF-2 can be used as a cell culture supplement to optimize the culture conditions for a broad variety of cell types, such as marrow stromal cells (MSCs), neural cells, and endothelial cells.

Optimal concentration for a specific application should be determined by a dose-response experiment.

2. Reference

1. Robinson, C. J. and Gaines-Das, R. (1994) The international standard for basic fibroblast growth factor (FGF-2); comparison of candidate preparations by in vitro bioassays and immunoassays. *Growth Factors* 11: 9–16.

All protocols and data sheets are available at www.miltenyibiotec.com.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

MACS is a registered trademark of Miltenyi Biotec GmbH.

Copyright © 2008 Miltenyi Biotec GmbH. All rights reserved.

