



Annexin V Binding Buffer (20× Stock Solution)

Order no 130-092-820

Index

1. Description

1.1 Background and product applications

2. Product dilution

1. Description

Product format	25 mL Annexin V Binding Buffer (20× Stock Solution).
Product size	For the preparation of 500 mL Annexin V Binding Buffer.
Storage	Store at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

1.1 Background and product applications

The Annexin V Binding Buffer (20× Stock Solution) is used for the preparation of a 1× Annexin V Binding Buffer. The Annexin V Binding Buffer has been specifically designed for use with all MACS® Annexin V products.

Product applications

- The Annexin V Binding Buffer (20× Stock Solution) is recommended for any labeling and washing step using fluorochrome- or biotin-conjugated MACS® Annexin V, the Annexin V MicroBead Kit (# 130-090-201), or the Dead Cell Removal Kit (# 130-090-101).

2. Product dilution

For preparation of 10 mL 1× Annexin V Binding Buffer, dilute 500 µL of the Annexin V Binding Buffer (20× Stock Solution) with 9.5 mL of sterile, distilled water.

Alternatively, prepare 1× Annexin V Binding Buffer by diluting 25 mL of the 20× Stock Solution with 475 mL of sterile, distilled water. Store at 2–8 °C.

▲ **Note:** Handle under sterile conditions.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

MACS is a registered trademark of Miltenyi Biotec GmbH.

© 2006 Miltenyi Biotec GmbH. Printed in Germany.

